



Shanghai Genomics, Inc.
No. 647 Songtao Road
Shanghai 201203, China
Phone: 86-21-50802786
FAX: 86-21-50802783
www.shanghaigenomics.com

PolyTag-GFP

CATALOG NUMBER: SG3000-04

LOT NUMBER: 30000911-03

Molecular weight: 74kDa

Package Size: 500µl

DESCRIPTION: Recombinant protein containing His, HA, Flag, GST, Myc tags and GFP, used for Western Blot

Source: Purified from E. coli expressed recombinant PolyTag-GFP

APPLICATION: Western Blot

Period of Validity: 5 years

Buffer: SDS loading buffer (50mM Tris pH6.8, 2% SDS, 10% Glycerol, 0.1% bromophenol blue, 1.5%DTT and 4M Urea)

Concentration: 10 ng/uL

Usage: Loading 5 uL/lane is recommended. For research purpose only.

STORAGE/Handling: Protein is recommended being stored at -80°C. Avoid repeated freeze-thaw cycles.

Amino Acids Sequence:



