



Shanghai Genomics, Inc.
No. 647 Songtao Road
Shanghai 201203, China
Phone: 86-21-50802786
FAX: 86-21-50802783
www.shanghai-genomics.com

Monoclonal Anti-HA-Peroxidase (HRP) conjugated

| | |
|---------------------------------------|---|
| CATALOG NUMBER: | SG4130-30 |
| QUANTITY: | 200 μ l |
| ISOTYPE:: | Mouse IgG2b |
| DESCRIPTION: | Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (YPYDVPDYA) coupled to KLH. |
| APPLICATION: | Western Blot |
| FORMAT: | 1 \times PBS with 50% glycerol |
| RECOMMENDED ANTIBODY DILUTION: | Western Blot: 1:500 |



Figure 1: Western Blot analysis of cell lysate (lane 1) and cell lysate expressing HA-tagged protein (lane 2).

| | |
|--------------------------|--|
| STORAGE/HANDLING: | Antibody is recommended being stored at -20 . Avoid repeated freeze-thaw cycles. |
| USAGE: | FOR RESEARCH USE ONLY. NOT FOR HUMAN USE. |

Background:

Epitope tags consisting of short sequences recognized by well-characterized monoclonal antibody have been widely used in study of protein expression in various systems. Our anti- HA-peroxidase conjugated antibody provides a simple solution to detect the expression of protein with a HA tag.